# ClinCal® - Calibrator

## Serum Calibrator Set lyophilised / Serum-Kalibrator-Set lyophilisiert

### FOR AMIODARONE / FÜR AMIODARON

#### Intended use:

ClinCal® serum calibrators are used for calibration of the analysis-system. These lyophilised calibrators are based on calf serum and are available as a calibrator set with 4 levels, including a blank. The calibrators have to be prepared like patient samples in one series of analyses.

#### Reconstitution:

Add exactly 2.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

#### Storage and stability:

Originally closed and stored at 4 °C in the dark, the lyophilised serum calibrators are stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 8 hours when stored at + 25 °C in the dark
- at least 5 days when stored at + 4 °C in the dark
- at least 3 months when stored at - 20 °C in the dark (avoid repeatedly freezing and thawing)

#### Notes:

The mean values of the analytes are chosen in a way that a linear calibration curve can be generated over a suitable range of concentration. According to quality assurance all ClinCal® calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 0.79 %.

#### Mean values:

The mean values were established in independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council (RiliBäk) with statistic methods.

#### Pack size:

ClinCal® Serum Calibrator Set for Amiodarone  
Level 0 - 3  
4 x 1 x 2 ml, **order no.: 9909**

#### Precautions:

This ClinCal® Calibrator Set was manufactured from pooled calfserum. It should be considered as potential infectious and treated with appropriate care.

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### Level 0

This blank-level is manufactured from the same pool of calf serum like the others, but it was not spiked with Amiodarone or Desethylamiodarone. Therefore the concentration of the analytes is zero.

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CHEMICALS + INSTRUMENTS GmbH      Labortechnik