ClinChek[®] - Control

Serum Control lyophilised / Kontrollserum lyophilisiert

FOR 25-OH VITAMINS D2/D3 / FÜR 25-OH VITAMINE D2/D3

Intended use:

ClinChek[®] serum controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on calf serum and are available with mean values in two different ranges of concentration. The controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 2.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved the solution is ready to use.

Storage and stability:

Originally closed and stored below -18 °C in the dark, the lyophilised serum control is stable for 36 months, but not longer than the expiration date printed on the label. After reconstitution the stability of the analytes is:

- at least 3 days when stored at 15 - 30 °C in the dark
 - at least 7 days when stored at 2 - 8 °C in the dark
 - at least 3 months when stored below -18 °C in the dark
- (avoid repeatedly freezing and thawing)

Notes:

The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek[®] controls have to pass strict qualitycontrol procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 1.81 %.

Mean values:

The mean values and confidence intervals were established in additional independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council (RiliBäk) with statistic methods.

Pack size:

ClinChek[®] Serum Control:
Level I
10 x 2 ml, **order no.: 35080**
Level II
10 x 2 ml, **order no.: 35081**
Level I, II
2 x 5 x 2 ml, **order no.: 35082**

Precautions:

The calf serum that was used for manufacturing the controls should be considered as potentially infectious and treated with appropriate care.

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CHEMICALS + INSTRUMENTS GmbH Labortechnik