ClinChek[®] - Control

Serum Control lyophilised / Kontrollserum lyophilisiert

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Intended use:

ClinChek[®] serum controls are used for internal quality assurance in the occupational and environmental toxicological laboratory as well as in the clinical-chemical one. These lyophilised controls are based on human serum and are available in two different ranges of concentration. After reconstitution, the controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 3.0 ml of analytical grade or bidistilled water to the vial and close it well. Then dissolve the lyophilised material completely. For this purpose we recommend a so-called roller mixer; a treatment of 30 min with such an instrument is sufficient. Otherwise incubate the samples for 1 to 2 hours at room temperature by shaking it carefully from time to time. Prior to use mix the controls well again. In case there are any precipitates, please centrifuge.

Storage and stability:

Originally closed and stored at 2 - 8 °C the lyophilised serum control is stable for 48 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 8 hours when stored at 15 - 30 °C
- at least 4 days when stored at 2 - 8 °C
- at least 1 month when stored below -18 °C (avoid repeatedly freezing and thawing)

Notes:

Routine analyses of the controls are performed according to the testprocedure for patient samples. The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek[®] controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 1.31 %.

Mean values:

The mean values and confidence intervals were established in independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council (RilBÄk) with statistic methods. For the determination of the trace elements, atomic absorption spectrometry (AAS) was used. Even if you use different techniques of analysis your results should be within the given ranges of the data sheet.

Pack size:

ClinChek[®] Serum Control
Level I
10 x 3 ml, **order no.: 8880**
Level II
10 x 3 ml, **order no.: 8881**
Level I and II
2 x 5 x 3 ml, **order no.: 8882**

Precautions:

The human serum that was used for manufacturing the controls was tested for the following infectious markers and found negative: HIV 1/2- and HCV-antibodies, hepatitis B-surface antigen. Nevertheless, the serum controls should be considered as potentially infectious and treated with appropriate care.

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CHEMICALS + INSTRUMENTS GmbH