## ClinChek® - Control

### Urine Control lyophilised / Kontrollurin lyophilisiert

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##### Intended use:

ClinChek® urine controls are used for internal quality assurance in the occupational and environmental toxicological laboratory as well as in the clinical-chemical one. These lyophilised controls are based on human urine and are available in two different ranges of concentration. After reconstitution, the controls have to be prepared like patient samples in one series of analyses.

##### Reconstitution:

Add exactly 5.0 ml of analytical grade or double distilled water to the vial and close it well. Then dissolve the lyophilised material completely. When all material is dissolved, the solution is ready to use.

##### Storage and stability:

Originally closed and stored at 2 - 8 °C the lyophilised urine control is stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 8 hours when stored at 15 - 30 °C
- at least 7 days when stored at 2 - 8 °C
- at least 1 month when stored below -18 °C (avoid repeated freezing and thawing)

##### Notes:

Routine analyses of the controls are performed according to the test procedure for patient samples. The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.5 %. The average residual moisture of this lot is 0.69 %.

##### Mean values:

The mean values and confidence intervals were established in independent reference laboratories in cooperation with the **INSTITUT UND POLIKLINIK FÜR ARBEITS-, SOZIAL- UND UMWELTMEDIZIN DER UNIVERSITÄT ERLANGEN - NÜRNBERG**. The calculation of the mean values and confidence intervals was done according to the guideline of the German Medical Council (RiliBäk) with statistic methods. For the determination of all parameters, chromatographic techniques were performed. Even if you use different techniques of analysis your results should be within the ranges indicated on the data sheet.

##### Pack size:

ClinChek® Urine Control  
Level I and II  
2 x 5 x 5 ml, **order no.: 8931**

##### Precautions:

These controls were manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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RECIPE                      CHEMICALS + INSTRUMENTS GmbH