ClinChek® - Control

Whole Blood Diabetes Control lyophilised / Diabetes-Vollblutkontrolle lyophilisiert

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Intended use:

ClinChek® whole blood controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human whole blood and are available with mean values in the normal as well as in the elevated range.

Reconstitution:

Add exactly 0.5 ml of HPLC water to the vial. To dissolve the material completely, we recommend a so-called roller mixer. A treatment of 15 minutes on such an instrument is sufficient.

If it is necessary to have a defined buffer-constitution of the control for further sample preparation (e.g. for pH-adjustment), you can also dissolve the lyophilised whole blood control directly in the appropriate buffer.

Procedure:

The reconstituted controls have to be pre-treated the same as patient samples and analyzed in accordance with instructions of the used analytic method.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised whole blood control is stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analyte is:

- at least 8 hours when stored at 15 - 30 °C
- at least 7 days when stored at 2 - 8 °C
- at least 30 days when stored below -18 °C (avoid repeatedly freezing and thawing)

Notes:

Routine analyses of the controls are performed according to the testprocedure for patient samples. The concentrations of the analyte are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® controls have to pass strict quality-control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 1.59 %.

Total hemoglobin concentration of the controls is 14.8 g/dl (for both levels).

Mean values:

The mean values and confidence intervals have been established in independent reference laboratories with supervision of an institution that is recommended by the guideline of the German Medical Council for quality assurance (RiliBäk). The Hb A1c values of the different methods are calibrated in accordance with the IFCC reference method. The metrological traceability is within the responsibility of the corresponding manufacturer of the Hb A1c-test.

Pack size:

ClinChek® Whole Blood Diabetes Control for Hemoglobin A1c Level I, II
2 x 5 x 0.5 ml, **order no.: 11092**

Precautions:

The human whole blood that was used for manufacturing the controls was tested for the following infectious markers and found negative: HIV 1 and 2, HCV and hepatitis B-surface antigen. Nevertheless, the whole blood controls should be considered as potentially infectious and treated with appropriate care.

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The mean values were derived from an external quality assessment in January 2010. For conversion of the values into the NGSP-unit, please use the following formula: **NGSP [%] = IFCC [mmol/mol] x 0.09148 + 2.152**

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In Germany Hb A1c is subject to the guideline of the German Medical Council for quality assurance. For assessment of the internal quality control, please also refer to table B 1a of the guideline.

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CHEMICALS + INSTRUMENTS GmbH Labortechnik